

Order Form

email: purchase@fireflyfriends.com



GottaGo

GottaGo is a unique portable toilet seat for children with disabilities. GottaGo is designed for use in the home and wherever you travel. This consistent approach is important when toilet training children with additional needs. GottaGo supports children in the squat position, which is the optimal position for toileting. GottaGo is available in two sizes, covering children aged approximately 2-9 years old.

Dealer		Marked for		Account	
Date		PO		Ordered by	
Phone		Email			
Ship to address					

Product Description	Product Codes	Product Chosen	MSRP
GottaGo Size 1 - Includes Seat, Frame, and all accessories	183-1600		\$1,895
GottaGo Size 2 - Includes Seat, Frame, and all accessories	183-2600		\$1,895
GottaGo Size 1 - Seat accessory only	183-616		\$699
GottaGo Size 2 - Seat accessory only	183-617		\$699

Seat	Size 1	Size 2
Age	2-5 years	4-9 years
Hip Width	8.8 inches	10.6 inches
Minimum Seat depth*	8.6 inches	11.4 inches
Backrest incline angle (forward)	8°	25°
Thigh angle (off horizontal)	40°	40°
Backrest height	9.4 inches	14.1 inches
Weight (seat and frame)	12.7 lbs	13.4 lbs

Frame	
Frame folded size	19.2 x 22.4 x 7.4 inches
Frame floor footprint	20.4 x 28.3 inches (at maximum height)
Frame floor footprint	15.7 inches
Maximum user weight	66lbs
Height (lowest part of the seat)	1 - 11.8in 2 - 18.1in 3 - 20.8in