



## WeGo

The lightweight push chair for kids with mobility issues

- Compatible with the Firefly GoTo Seat\*
- Lightweight
- Recline angle up to 140 degrees
- Full accessory pack including canopy and rain cover
- Wheels suitable for rough terrain
- Easy to fold
- Small footprint so can fit in most car trunks.

\*Please check sizing guide for exact measurements\*

Dealer		Marked for		Account	
Date		PO		Ordered by	
Phone		Email			
Ship to address					

Product Description	Product Codes	Product Chosen	MSRP
WeGo Push Chair	182-500		\$1,465

## WeGo Size Guide

Footprint Length (Collapsed)	46 inches
Footprint Height (Collapsed)	15 inches
WeGo Height	48 inches
Seat Depth	10 inches
Seat Height	20 inches
Seat Width	14 inches
Footrest Shortest Setting	10 inches
Footrest Longest Setting	13 inches
WeGo Weight	24 lbs
WeGo Weight with GoTo	29 lbs

ONE SIZE

MAX 66 Lbs

\*Please check sizing guide for exact measurements\*

