

# Order Form

email: [purchase@fireflyfriends.com](mailto:purchase@fireflyfriends.com)



## Scooot

**Scooot** is a configurable 4-in-1 mobility rider designed to help children with special needs play, explore and participate in a way that best suits their physical and cognitive abilities.

Dealer		Marked for		Account	
Date		PO		Ordered by	
Phone		Email			
Ship to address					

Product Description	Product Code	Product Chosen	MSRP
Scooot 4 in 1	164-620		\$1,575
Scooot 4 in 1 with Advanced Backrest	164-621		\$1,845

## Scooot Size Guide

Size	Approx age	Maximum user weight	Seat base width	Seat base length (from back to front of seat)
4 in 1	2-6 years	48.5 lbs	9 inches	15.5 inches

  

Length from hip to end of footplate when seated	Backrest width	Backrest height	Wheel arch width	Height of seat base off floor
24.5 inches	10.5 inches	12 inches	9 inches	6 inches

  

Laterals length	Lateral length settings	Measurements
Advanced backrest cushion accessory	16.5 inches	31 inches (min) 40 inches (max)
	Pooosh Pole	H 35 inches (max from floor) H 31 inches (min from floor)